

**Pretreatment and Determination of Chlortetracycline,  
Oxytetracycline, and Tetracycline in Chicken and Pork by  
SepLine/HPLC-UV**

LabTech, Inc.

Tetracycline antibiotics are used to prevent and cure animal diseases. However, long term usage of tetracyclines may cause residue accumulation in animal tissues, which affects the consumer's health. Both European Union and China have restrict rules about the amount limits of tetracyclines in animal tissues. Since fast-grown chicken, feed with certain drugs attracted great attentions recently, determination of drug residues from animal tissues is critical at this point. This present study describes the analytical methodology by using SepLine-1 automated solid phase extraction (SPE) system and high-performance liquid chromatography (HPLC) coupled with ultraviolet detector for the sensitive determination of antibiotics chlortetracycline, oxytetracycline, and tetracycline in chicken and pork. The method is simple-to-use and has high recovery, which can be widely used as routine tetracyclines testing.

## 1. Experimental

### 1.1 Instrumentation and Materials

SepLine -1 Automated SPE System (LabTech)

HPLC LC600 System (LabTech)

SPE Cartridge: LabTech PLS 150 mg 6 mL

Chlortetracycline, Oxytetracycline, and Tetracycline standards, purity  $\geq 95\%$

Sodium Hydrogen Phosphate Solution (0.2 mol/L): Dissolve 28.41 g sodium hydrogen phosphate in 1000 mL water.

Citric Acid Solution (0.1 mol/L): Dissolve 21.01 g citric acid in 1000 mL water.

Mcllvaine Buffer: Mixture of 1000 mL 0.1 mol/L Citric Acid and 625 mL 0.2 mol/L Sodium Hydrogen Phosphate Solution.

Na<sub>2</sub>EDTA-Mcllvaine Buffer (0.1 mol/L): Dissolve 60.5 g Ethylenediaminetetraacetic acid disodium salt (Na<sub>2</sub>EDTA) in 1625 mL Mcllvaine Buffer, and shake well.

Oxalic Acid Solution (0.01 mol/L): Dissolve 0.882 g oxalic acid in 700 mL water.

Mobile Phase: Acetonitrile + Methanol + 0.01 mol/L oxalic acid solution. (20%/10%/70%)

Chlortetracycline, Oxytetracycline, and Tetracycline Stock Mixed Standard Solution: Weigh 10 mg of each standards and dissolve into 100 mL Methanol.

Chlortetracycline, Oxytetracycline, and Tetracycline Mixed Working Standard Solution: Dilute Stock Mixed Standard Solution into different concentration (0.25 mg /L, 0.5 mg /L, 1.0 mg /L, 5.0 mg/L) by mobile phase.

Methanol (AR.)

**LabTech, Inc.**

Address: 114 South Street, Hopkinton,  
MA 01748, U.S.A

Tel: (508) 435-5500

Fax: (508) 435-5595

Oxalic acid (AR.)

Pure Water

## 1.2 Sample Preparation

### Sample Extraction

- Mince and mix animal tissues, then transfer into a clean container. Weigh 6 g samples to the nearest 0.01 g, then transfer to 50 mL centrifugal tubes with 30 mL 0.1mol/L Na<sub>2</sub>EDTA-McIlvaine buffer. Mix sample for 1 min on shaker and then vortex for 10 min.
- Centrifuge for 10 min at 5000 rpm. Transfer the supernatant to another centrifugal tube with additional 20 mL buffer. Centrifuge again and combine the supernatants.

### Method Summary

Table 1. Extract method of Tetracyclines determination in animal tissues by SepLine-1 Automated SPE system. (SPE Cartridge: LabTech PLS 150 mg 6 mL)

Step	Solvent	Volume (mL)	Flow Rate (mL/min)	Time (Sec)
Prewet 1	Water	5.0	5.0	0
Prewet 2	Buffer	5.0	5.0	0
N <sub>2</sub> Purging	NULL	0.0	0.0	10
Sample Loading	Sample	10.0	2.0	0
N <sub>2</sub> Purging	NULL	0.0	0.0	180
Flow Path Rinse 1	Water	5.0	5.0	0
N <sub>2</sub> Purging	NULL	0.0	0.0	180
Elution	Methanol	10.0	2.0	0
Elution	Methanol	10.0	2.0	0
N <sub>2</sub> Purging	NULL	0.0	0.0	100

The eluate was dried by rotatory evaporator at 40 °C, and then reconstructed into 1 mL mobile phase.

## 1.3 Apparatus

### HPLC Method

Column: LabTech C18 (250 mm × 4.6 mm, 5µm)

Column Temperature: 30°C

Mobile Phase: Acetonitrile + Methanol + 0.01 mol/L Oxalic Acid (v/v/v = 20%/10%/70%)

Flow Rate: 0.5 mL/min

Sample Volume: 20µL

Wavelength: 350 nm

**LabTech, Inc.**

Address: 114 South Street, Hopkinton,  
MA 01748, U.S.A

Tel: (508) 435-5500

Fax: (508) 435-5595




Figure 1. Chromatogram of 1 ppm Chlortetracycline, Oxytetracycline, and Tetracycline Standards in Chicken after SPE.

## 2. Results and Discussion

**LabTech, Inc.**

Address: 114 South Street, Hopkinton,  
MA 01748, U.S.A

Tel: (508) 435-5500

Fax: (508) 435-5595






Figure 2. Calibration Curves of Chlortetracycline, Oxytetracycline, and Tetracycline.

**LabTech, Inc.**

Address: 114 South Street, Hopkinton,  
MA 01748, U.S.A

Tel: (508) 435-5500

Fax: (508) 435-5595

The linear range is between 0.25 mg/L – 5.0 mg/L. The equations are listed below, where x is sample concentration, y is peak area.

Oxytetracycline	$y = 80317x + 285.7$	$r^2 = 0.999$
Tetracycline	$y = 60467x + 9499$	$r^2 = 0.999$
Chlortetracycline	$y = 30913x + 3524$	$r^2 = 0.999$

9 mL animal tissue samples were spiked with 1 mL 1 ppm of standards. Then these spiked samples were purified through SepLine Automated SPE system, and then analyzed by HPLC. The mean recoveries of chlortetracycline, oxytetracycline, and tetracycline were listed in table 1 and table 2.




Figure 3. The LC chromatogram of chicken samples spiked with 1 ppm after SPE.

Table 1. Recovery of Spiked Chicken Samples.

	Compound	Sample 1		Sample 2		Sample 3	
		Measured (ppm)	Recovery (%)	Measured (ppm)	Recovery (%)	Measured (ppm)	Recovery (100%)
1	Oxytetracycline	0.874	87.4	0.882	88.2	0.979	97.9
2	Tetracycline	0.90	90.0	0.892	89.2	0.927	92.7
3	Chlortetracycline	0.994	99.4	0.897	89.7	0.971	97.1

**LabTech, Inc.**

Address: 114 South Street, Hopkinton,  
MA 01748, U.S.A

Tel: (508) 435-5500

Fax: (508) 435-5595

Table 2. Recovery of Spiked Pork Samples.

	Compound	Sample 1		Sample 2		Sample 3	
		Measured (ppm)	Recovery (%)	Measured (ppm)	Recovery (%)	Measured (ppm)	Recovery (100%)
1	Oxytetracycline	0.839	83.9	0.963	96.3	0.979	97.9
2	Tetracycline	1.016	101.6	1.08	108	0.906	90.6
3	Chlortetracycline	0.981	98.1	1.09	109	1.09	109




Figure 4. The LC chromatogram of chicken sample after SPE system. (Tetracyclines are not detectable)

### 3. Conclusion

This study develops the method by using SepLine Automated SPE system and HPLC to determine chlortetracycline, oxytetracycline, and tetracycline in chicken and pork. The results demonstrated that this method is fast, reliable and precise for chlortetracycline, oxytetracycline, and tetracycline determination in animal tissues.

**LabTech, Inc.**

Address: 114 South Street, Hopkinton,  
MA 01748, U.S.A

Tel: (508) 435-5500

Fax: (508) 435-5595